



iMWard EPD Bedside / Door Sign Patient Card


Real-time reminder to reduce workload and improve care quality



- Integrate patient information from hospital information system and display it on the iMWard Bedside Patient Card automatically.
- Provide real-time medical care reminders to improve care quality and interactions between caregivers and patients.
- Reduce workload on paperwork - first step towards paperless hospital.





The advantage of iMWARD Bedside Patient Card



- ✓ Display pivotal information. Optimize communications between caregivers and patients.
- ✓ Integrate with hospital information system. Update patient information instantly.
- ✓ Reduce caregivers' workload on paperwork.
- ✓ Achieve green hospital with paperless iMWARD Bedside Patient Card solution.
- ✓ Patient-care data analysis.


*No light pollution
and low power consumption*

Optimize workflow with iMWARD Bedside Patient Card



Item	Traditional bedside card	iMWARD Bedside Patient Card
Presented information	Limited information	Customized information
Information update	Manually	Instantly
Human error	High	No
Emergency call	No	Yes

iM Ward Bedside Patient Card System Architecture



iM Ward Patient Management Platform

The most convenient way - manage your epaper in one place

01 Patient management

personalized information such as treatment process

02 Hospitalization management


Notification query, announcement, account authority

03 Device monitoring

battery status, connection status and connection time



Product Sizes – From 9.7" to 42"



9.7"	13.3"	13.3"	42"
Black/ White/ Red	Black/ White/ Red	60k color	Black/ White
4 physical buttons	4 physical buttons	-	capacitive touch
built-in 3100mAh	built-in 3100mAh	built-in 4900mAh	Output: 60W (12V/5A)
WiFi/BT/NFC	WiFi/BT/NFC	WiFi/BT	Lan
9.8"x7.9"x0.7"	12.4"x10.2"x0.8"	11.8"x9.1"x0.9"	36.2"x27.6"x1.8"
LED indicator	1.HIPAA compliant 2.LED indicator 3.Battery swappable	-	support pain management

Why choose imedtac?

We are Taiwan's Healthcare AIoT expert, who can help you enable your medical AIoT business



1000+
bedside cards



10+
institutions

